


for deep uterine horn embryo transfer...

Cow Embryo Transfer Catheter *mo-No.4*




- Deep uterine embryo transfer catheter
- For 0.25mL embryo transfer straw
- Packing : Individually packed
- Package : 20 sets/box




「mo-No.4」 ensures successful embryo transfer to the deep part of uterus, is realized in cooperation with ZEN-NOH Embryo Transfer Center Japan.

- ◎ Insert “mo-No.4” into uterine cornu and plunge the connector. Flexible plastic tube will extend into the deep part of the uterus to release the embryo.
- ◎ Cover tube to prevent dirt is include in the set, and “mo-No.4” is individually packed. It can be taken out together and ready for single use.



How to use 「mo-No.4」

- [1] Open the individual pack from the connector side and insert the straw into the hole in the connector.
- [2] With a thin stick, insert the embryo and fluid from the straw into the tube inside the catheter. (When doing this it decreases the resistance if you cut the plug into half which makes it easier to push.)
- [3] Keep the catheter horizontal in order to prevent the contents from spilling, pull off the straw and connect a syringe (1mL or 2 mL) to the connector. (Syringe plunger should be pulled at this point to the amount to release the content. **ATTENTION : If the syringe is disconnected the content spills out.**)
- [4] With cover tube still on, take out the catheter from the packing bag and operate on the transfer.
- [5] When the head tip is inserted into the uterine cornu, push in the connector to extend the head tip. When it has reached the desired depth, push the syringe plunger slowly (taking 4~5 seconds) to release the embryo.



※Caution: Store in a cool and dry place, away from high temperature, humidity and direct sunlight. When the packet is wet or broken, it shall not be used. **This item is designed for single use only.**

MISAWAMEDICAL INDUSTRY CO., LTD.

351, Asahi-machi Kasama-shi, Ibaraki-ken, 309-1717 Japan

Phone: +081-296-77-8804 Fax: +81-296-77-8849

e-mail: sales@misawa-medical.co.jp <http://www.misawa-medical.com>