

Ggauge	OD	length											
		4mm	6mm	8mm	12mm	16mm	25mm	30mm	40mm	50mm	75mm	100mm	125mm
15G	1.80mm						●						
16G	1.60mm						●						
17G	1.40mm							●					
18G	1.20mm						●			●	●	●	●
19G	1.10mm												
20G	0.90mm									●			
21G	0.80mm						●			●			
22G	0.70mm							●					
23G	0.60mm						●			●			
24G	0.55mm						●						
25G	0.45mm						●			●			
26G	0.40mm				●	●		●		●			
27G	0.40mm	●	●	●	●	●	●			●			
30G	0.30mm	●	●	●	●		●		●				
31G	0.25mm	●	●	●	●		●						
32G	0.23mm	●	●	●	●								
33G	0.20mm	●	●	●	○								
34G	0.18mm	○	○	○									

* We are able to handle anything from supply of our standard size injection needles on an OEM basis to customer requirements for specified dimensions.

* Manufacture of special sizes of injection needles is accommodated through sample manufacture.

● STANDARD SIZES
○ UNDER REQUEST SIZES